
# MUHAN Co., Ltd. COMPANY PROFILE




## MANAGEMENTIDEOLOGY

To become a global company with the goal of growing into a global company, we will work together in every way to become an original, creative, and creative company, and grow into an international level manufacturing and distribution company.



· Employees' loyalty and ownership

## ORGANIZATIO



## COMPANY HISTORY



### History (2012~now)

- 12. 2012 Establishment of MUHAN Corporation, Inc.
- 12. 2012 AUDI, LENOMA Duty Free Shop and Department Store
  02. 2018 Henry Cotton License Agreement open a store
- 12. 2012 Nixon Online Distributor Agreement
- 12. 2012 Kenneth Cole National Distributor Agreement
- 11. 2015 Acquisition of the FOCE trademark
- 04. 2016 Swiss 150-year-old Brand Henry Blanc Acquire Global Trademark Rights
- 09. 2016 SCORPION trademark acquisition

- 12. 2017 License Tottenham Hotspur Collaboration Acquisition
- 02. 2018 CLAREO OEM Manufacturing Contracts
- 05. 2018 CODEMUNTZ OEM manufacturing contract concluded
- 09. 2018 Participation in the Hong Kong Trade Fair
- 01. 2019 MPGUARD protective clothing launching
- · 06. 2019 Establish a distribution department for Ili Coffee
- 09. 2019 Participation in the Hong Kong Trade Fair
- 03. 2020 Defense Clothes/FDA Approved

## PRODUCT DESCRIPTIONS



### medical protective clothing

medical gown

## PRODUCTION LINE



## PRODUCTION LINE







Protective Medical Cloth - Level D MP60G Related product proposal





Artificial blood permeability resistance

STEP 1~6 PASS

liquid penetration test(%)

MD penetration index 0.0% CD penetration index 0.0% ISO16603 Certification

## **MP GUARD is**


**«ΚΔΤ**Π

Level D protective clothing that effectively protects against the penetration of viruses, bacteria, blood and foreign substances from outside. The non-woven laminated fabric, which maintains moisture and waterproofing while delaying or blocking the penetration of various types of viruses, makes it more pleasant and durable for the caregiver (medical staff).



## **KATRi TEST PASS**

Completed blood and body penetration resistance test using synthetic blood Same as CE certified EN14126 / ISO 16603 certified




시 혐 항 목	시 험 겸 과	
~ 8 8 7	人員1	
인공혈액침투저항성 : KS K ISO 16603:2012		
Step 1	Pass	
Step 2	Pass	
Step 3	Pass	
Step 4	Pass	
Step 5	Pass	
Step 6	Pass	

#### ▲ Results of blood tests using synthetic blood

11 21 21 21	시 혐 결 3
시 혐 향 목	서료1
액체 침투시험 (%) : 고용노동부고시 제2020-35호, 별표 8의4	
u	
MO	
- 吾年지수(%)	0.0
- CLASS	3
00	
	0.0
- 침투지수(%)	3

▲ Results of Liquid Penetration Resistance Test

# Type of Protective clothing



# Product Features

Specialized SF fabric blocks nasal and

blood, and P.E coating makes you feel more secure.

04

SF Fabric / P.E Coating

02



#### Safe Boots cover

First, block off contaminants with banding treatment. It's safer because you can tie it with a string.



## Safety with double-sealed seal

primary blocking with zipper. Use the adhesive tape on the top of the zipper the second time. a safer protective suit.

# Banding to prevent

intrusive interception

Banding the wrists and ankles to prevent internal penetration of contaminants.



# Chart for protective clothing size




Size	Ť	Ť	Ť	İ
	М	L	XL	XXL
Available height	150-165	165-175	175-190	190-200
Total length	198	200	204	208

## About C/T

Product Type	Size Type	Box size	Weight	Quantity
MP60G	M/L/XL/XXL	570x400x630	16Kg	50sets(1box)
Product Type	Size Type	Box size	Weight	Quantity
Boots cover	FREE	570x400x630	13Kg	300sets(1box)

## Korea Clothing Testing Laboratory (KATRi) Test report



# Thank you

MUHAN Co., Ltd. CEO Lim, Sung-il

B1, 36-10, Sangdo-ro 34gil, Dongjak-gu, Seoul, Republic of Korea Mobile. +82-10-5064-7748 E-mail. minipos2@nate.com Tel. +82-2-3663-7377 Fax. +82-2-3663-8377